

## Product Data Sheet

### Prestige Plus – Machine Detergent Capsules

Chlorinated powder. Super concentrated capsule style detergent for use through dishwasher. Ideal for removing heavy food soils.

---

#### Product Description

Prestige Plus is a white chlorinated detergent, specifically designed for the use in your Ecowash automatic 'powder feed' detergent system. It is formulated to remove heavy food soils and tannin in all water conditions. Top of the line capsule technology. Super concentrated with high levels of chlorine to achieve the best results at lowest usages.

#### Characteristics

- White powder capsules
- Corrosive
- Non Flammable
- pH 14

#### Tools / Safety Equipment

- \* Goggle
- \* Rubber or PVC gloves.

#### Safety

This is a corrosive product and may cause severe burns. Risk of serious damage to eyes. Avoid contact with skin and eyes. Wear suitable gloves and eye/face protection when handling this product.

#### Application Procedure

1. Read the Safety Directions and Directions for Use, as well as Material Safety Data Sheet on Prestige Plus. This system should

only be used by trained responsible staff.

2. Remove empty capsule and dispose of appropriately. Do not reuse or refill capsule with other powders as this may pressurise and become extremely hazardous.
3. Lift the tab to expose screen on (label will be upright) capsule, do not unscrew lid or damage the screen inside lid. Insert capsule into dispenser bowl, ensuring label is turned to the front of the dispenser to enable reading of same.
4. Prestige plus will work at its best at a temperature of 50-60° C. do not use if water temperature into capsule exceeds 60° C as product may pressurise and become extremely hazardous.
5. Ensure Prestige Plus capsule has product in it at all time.

#### First Aid

- If poisoning occurs, contact a doctor or Poisons Information Centre (Ph: 131 126)
- If swallowed, do NOT induce vomiting. Give a glass of water
- If skin contact occurs, remove contaminated clothing and wash skin thoroughly.
- If in eyes, hold eyes open, flood with water for at least 15mins and see a doctor.

#### Pack Sizes

3 x 5kgs